HUR Cardio Line 800





HUR Cardio Line 800 med

Based on time-proven components such as e. g. the grid independent, highly precise IBS brake system, we have developed a new product line, which meets the highest standards in terms of quality and design. Together with Emotion Fitness HUR created a cardioline equipped with high class technology SmartCard.

SmartCard - Automated Workout

All machines of the HUR cardio line 800 med are equipped with a SmartCard reader as a standard feature. SmartCard automates the entire workout - just insert the card and start exercising.

Versatile Feedbacks

Via SmartCard, individual information for a fitness test or the workout data can be saved and transported to a SmartCard Software. After exercise is done, it gives versatile progress reports which enhance the user's motivation.

Design

HUR cardio line 800 med offers the continuous, strongly dimensioned and premium anodised aluminium profile. The plastic covers were reduced to the minimum and the LED-background light in the display ensures a higher comfort when reading the data and operating the machines.

Oy HUR Ab
Health & Fitness Equipment
Patamäentie 4 • FI 67100 Kokkola
Tel. +358 (06) 832 5500
Fax +358 (06) 832 5555
www.hur.fi





HUR sprint 800 SL/SE med

with the well known and proven lamella technology ensures supreme cushioning, smooth synchronization and a very low power consumption.

On the versions with inclination, the pulse controlled programs can be adjusted either via the variation of speed or inclination.

FUNCTIONS

- · Easy to handle, backlit display
- Very comfortable running feeling and optimal cushioning
- · Coded Polar® pulse receiver · Speed: 0–17 km in 0.1 km/h
- · Speed: 0–17 km in 0.1 km/h steps
- · 0-15% incline
- . Heart rate controlled programs either incline or pulse controlled
- Up to 1 km: 1 m precise odometer, then 10 meter precise



HUR relax 800 - The recumbent ergometer ensures a workout in a comfortable position and is therefore perfectly suited for long training sessions as well as heavy users. For those, the free step through is mandatory.

For an even more comfortable horizontal adjustment, with this machine the drive unit rather than the seat is being moved.

FUNCTIONS

- Easy to handle, backlit display
- · Easy step through
- Gas spring supported, horizontal adjustment of drive unit
- Grid independent
- · Coded Polar® pulse receiver



HUR body 800 - The upper body ergometer can be easily adjusted via a gas spring support. It can be adjusted for the seated or standing training. The platform is only 5 mm high and can therefore easily be accessed by wheelchairs, the seat can be removed completely.

The optional length adjustable crank arms can be put into a synchronous mode in order to simulate a handibike use.

FUNCTIONS

- · Easy to handle, backlit display
- Forward and backward movement possible
- · Gas spring supported, stepless height adjustment of the drive unit for ideal positioning
- · Horizontal seat adjustment
- · Seat removable for the standing workout, use of Pezziball or for wheelchair access
- · Grid independent
- · Coded Polar® pulse receiver



The HUR cross 800 med

impresses with its upright training position, is demanding for a large number of muscle groups and is gentle to the joints. The elliptical movement ensures a smooth movement with a good balance of arm and leg movement.

At the end or interruption of every training, the movement will be stopped completely, before the Cool Down starts.

FUNCTIONS

- · Easy to handle, backlit display
- Grid independent
- Coded Polar® pulse receiver
- Forward and backward movement possible
- · Big, anti-slip footsteps
- Rpm limit for extra safety



The HUR cycle 800 is the classic machine of cardio training and the ideal test and training ergometer for the professional user. The grid independence of the unit can be assured even with most of the various options, which guarantees flexibility and saves money.

For the health oriented use, the HUR cycle 800 offers an easy step through as well as an ergonomic training position.

The optional length adjustable crank arm is an extremely important feature for all users with limitations in the range of motion.

FUNCTIONS

- Easy to handle, backlit display
- · Easy step through
- · Horizontally and vertically adjustable saddle
- · Grid independent
- Coded Polar® pulse receiver







HUR line 800 med General Information

- Optional hand pulse system
 (except for the HUR body 800 med)
- Length-adjustable crank arm available
- Power independent, 25-500 watts in rpm independent mode, 25-1000 watts in rpm dependent mode (except for HUR sprint 800 SE/SL med)
- Medically approved by the CE-Label for medical products EN 93/42 EEC
- Compatible with HUR SmartCard system

For more information, please contact us at www.hur.fi



HUR cycle 800 med

- * Easy step through
- * Various interfaces to be adapted to different external systems
- * To be used as a test ergometer as well
- * Size (LxWxH): 130x65x146 cm
- * Weight: 74 kg



HUR relax 800 med

- * Easy step through
- * Horizontal adjustment of the drive unit
- * Optional user weight of 250 kg
- * Size (LxWxH): 159x65x137 cm
- * Weight: 82 kg



HUR sprint 800 SL/SE med

- * 0-17 km/h in 0,1 km/h steps, 0 to 15% incline (HUR sprint 800 SE med without incline)
- * Energy saving, brushless 1100W DC Motor, toothed belt lamella system
- * Integrated Emergency-Stop-System
- * Options: increased max speed, reversed direction of movement
- * EN60601-1: 230V/50-60Hz/EWG, separate fuse protection required
- * Size (LxWxH): 198x79x145 cm
- * Weight: 170 kg



HUR body 800 med

- * Forward and backward movement
- * Easy, stepless height adjust ment for the seated or standing workout
- * Seat easy to remove
- With length adjustable crank arm also handibike simulation possible
- * Size (LxWxH): 122x80x174 cm
- * Weight: 108 kg



HUR cross 800 med

- * Forward or backward movement
- * All body movement
- * Compensation of the seated worklife of many users
- * High quality slide and ball bearings for long product life
- * Size (LxWxH): 215x65x157 cm
- * Weight: 118 kg

